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How to use your Kinepod active insoles ?

There are two ways to use your Kinepodic insoles based on the chosen action :

1. Stimulating action :

Facilitating your muscles' work when walking and reducing joint constraints. Wear your kinepodic insoles when walking to encourage a physiologically plantar movement distributed across your entire body. Your insoles will enable you to train with the minimum desirable mobility for your body. Thanks to your insoles your body's muscular tensions will even out (the muscles that are overworking will reduce their activity and the muscles doing little work will increase their activity) and the impact coming from the ground will be better absorbed. To achieve a muscular balancing effect, it is recommendable to wear Kinepod insoles at least 1 hour per day. To achieve an amortisation effect, it is recommendable to wear Kinepod insoles whenever you walk or run. Respect the recommended adaptation time to allowing weak muscles to strengthen without risking injury. Allow your body the time it needs to balance itself.

2. Postural Action :

Relieve muscle and joint tensions using a few simple exercises on your Kinepod orthopaedic insoles each morning and during the day. Remove your kinepodic insoles from your shoes and place them on a hard floor (avoid carpets). Stand on the insoles. You are now ready for the simple exercises.

- Turn your head to each side, left and then right, and bow it on each side.
 - Turn your torso (including shoulders) to each side and bow it on each side.
 - Keeping your legs straight, lean forward as if you wanted to touch the floor with your hands.
 - Crouch by bending at the knees then straighten back up.
- Repeat the exercise two or three times as you wish. Your tensions should be relieved immediately.

Side effects

In the case of pains under the foot, it is preferable to reduce the time during which you wear the insoles. It is likely that the bearing has been too intense. If your pain persists, consult your doctor or a specialist.

Contraindications

Bearing the insoles is contraindicated with anti-pronation and anti-supination shoes, as well as with shoe having reinforcements for the arch of the foot, sandals or open shoes.

Tips for your insoles

1. Wear suitable shoes to achieve the maximum benefit from your soles It is recommended to wear neutral shoes (so-called universal shoes) for sport.

2. Shoes

The ideal shoe is: Light with a flat external sole (heelless and 0.5 to 1 cm thick) ; Closed at the instep (velcro or laces) As broad as your foot (at the heel, midfoot and forefoot) ; Semi-rigid under the heel (hard rubber) ; Supple under the forefoot ; With a hard buttress behind the heel. At home, it is recommended to be barefoot or wear only fine slippers, closed at the rear. Avoid shoes open at the back (such as flip-flops). These shoes are not designed for lengthy walks.

Sports shoes: it is essential to use neutral sports shoes (so-called universal shoes) to ensure the best performance of your sole. Given the variation of the shoes quality, Kinepod cannot be taken responsible in any case in which the shoes might be damaged by the insole.

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Given the variation of the shoes quality, Kinepod cannot be taken responsible in any case in which the shoes might be damaged by the insole.

3. Instructions to insert your soles into your shoes
If your shoes contain removable internal soles: Remove the removable sole before inserting the Kinepod orthopaedic sole in its place. If the Kinépod orthopaedic sole is too big, place the original sole on top of the Kinepod orthopaedic soil. Use a marker to trace the shape of the removable sole and trim the surplus from the new sole where necessary.

If your shoes do not have removable soles: it is recommended to insert the Kinepod orthopaedic sole in the shoe and to check whether it curls up at the front. If so, remove the sole and trim 2-3mm from the front in a round shape. Repeat until the sole fits flat inside the shoe.

4. Allow time to adapt to wearing your orthopaedic soles Wearing the sole requires an adaptation period which varies from one person to another. Generally it is between 3 and 8 weeks. Wear the soles progressively when walking for the 2 first weeks, i.e. 1 to 2 hours the first day, 2 to 3 hours the second, and so until building up to the 15th day. The goal here is to eventually wear them continuously. You should only wear the soles during sport after the 2 weeks of adapting to walking. Nevertheless, you will need time to adapt to running, i.e. 5 to 10 minutes the first day, 10 to 15 minutes the second and so until building up to the 15th day.

5. Care
The soles can be washed in soapy water. However, do not place them on a heat source (also avoid washing machines). They should be dried with a towel and at room temperature. If your soles are slipping inside your shoes, it is recommended you use double-side, sticky-back velcro (available in all sewing shops). Normal wear for Kinepod soles: over time the hard base will rub through the sole's underside (velvet). This physiological wear in no way changes the soles' effect.

6. Things to avoid
Do not wear shoes that are unsuitable for your feet. Doing so may counteract the sole's therapeutic effect and be harmful to your feet. Do not attempt to speed up the adaptation period even if your soles seem comfortable. Do not try to modify a sole yourself (except trimming). In addition to the risks, doing so invalidates the sole's guarantee of origin. Do not use greasy substances on your soles. Never expose your soles to intense heat sources and avoid all contact with abrasives.

7. Warranty
- 6 months against breakage, crackling or fragmentation of the internal rigid central base.
- 6 months against material or labour defects. This warranty is rendered null and void if the sole has been subjected to non-compliant, unusual or abusive treatment, to intentional damage or modifications performed by a third party other than the KINEPOD laboratory.

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